

ACTUATED MEDICAL® FACE SHIELD

Protecting all + Guarding futures

Reusable facial covering for most environments.

Face coverings, such as face shields, help minimize respiratory droplet contact between people. Compared to face masks, face shields enable people to see one another's facial expressions while minimizing accidental face touching.


With multiple available models, all Actuated Medical face shield models use the same 0.02 inch thick clear polycarbonate. Additionally, reinforced headband and hard hat models allow room for glasses or goggles as well as masks.



Actual Face Shield may appear different than pictured.

Key Features

- + Meets ANSI Z87 standard
- + Clear polycarbonate for visibility
- + 300% thicker than typical face shields for reuse and durability
- + Open design reduces fogging
- + Chemical resistance to 1:8 bleach/water solution for over an hour
- + Easy to clean
- + Headbands may include Acrylonitrile Butadiene Styrene (ABS) and elastic



Standard



Small



Hard Hat

These face shields have not been FDA cleared or approved. The face shields have been authorized by FDA under an EUA for use by healthcare providers as personal protective equipment, and are only authorized for the duration of the declaration that circumstances justifying the authorization of emergency use under Section 564(b)(1) of the Act, 21 USC 360bbb-3(b)(1) unless the authorization is terminated or revoked sooner.

THE SHIELDS ARE PROVIDED "AS-IS," AND ACTUATED MEDICAL DISCLAIMS ALL WARRANTIES, INCLUDING ANY (A) WARRANTY OF MERCHANTABILITY; OR (B) WARRANTY OF FITNESS FOR A PARTICULAR PURPOSE; WHETHER EXPRESS OR IMPLIED BY LAW, COURSE OF DEALING, COURSE OF PERFORMANCE, USAGE OF TRADE OR OTHERWISE. ACTUATED MEDICAL SPECIFICALLY DISCLAIMS THAT THE SHIELDS WILL PREVENT INFECTION OR THE TRANSMISSION OF VIRUSES OR DISEASES, INCLUDING COVID-19/CORONAVIRUS. ACTUATED MEDICAL FURTHER DISCLAIMS THAT THE SHIELDS WILL PREVENT PERSONAL INJURY OR DEATH. BUYER BEARS THE ENTIRE RISK AS TO THE PERFORMANCE, OPERATION, AND QUALITY OF THE SHIELDS.

Buyer acknowledges and understands that the shields have not been industry tested and may cause injury. Buyer agrees to indemnify, defend, and hold Actuated Medical harmless against any and all liability or loss arising in conjunction with or resulting from use of the shields (including any and all attorneys' fees and expenses incurred by Actuated Medical).

ORDER ONLINE:

www.FaceShield.ActuatedMedical.com

For larger quantities:

Call +1 (814) 355-0003 ext. 117

Email FaceShield@ActuatedMedical.com



Innovative motion + Positive outcomes



WOMEN OWNED

1100793368-009
Pat. ActuatedMedical.com/ip